


PRODUCT SPECIFICATIONS

PRODUCT COMPONENTS



OPERATION PROCEDURE

Nasopharyngeal swab

Oropharyngeal swab

Anterior nasal swab

Unscrew the lid of the extraction tube

Rotate about 10 times, leave the swab in the extraction buffer for 1 minute

Squeeze the swab tip

Cover the lid

Unscrew the small cap at the top of the extraction tube

Add 80 µL (about 3~4 drops) to the sample well

15 mins

2019-nCoV Antigen Result Interpretation

Positive: C (red line), T (red line)

Negative: C (red line), T (no red line)

Invalid: C (no red line), T (red line)

PERFORMANCE

Reagents	PCR		Total
	Positive	Negative	
W196-2019-nCoV Antigen Test (Lateral Flow Method)	Positive	208	209
	Negative	4	365
Total	212	362	574

Sensitivity:
98.11% (95%CI: 95.24%~99.48%)

Specificity:
99.72% (95%CI: 98.47%~99.99%)

Total agreement:
99.13% (95%CI: 97.98%~99.72%)

Data for clinical evaluation*				Data for asymptomatic patient							
Reagent	Confirmed by Referenced method PCR			Total	Reagent	Confirmed by Referenced method PCR			Total		
	Positive	Negative				Positive	Negative				
W634-2019-nCoV Antigen Test (Lateral Flow Method)	Positive	66	5	71	W634-2019-nCoV Antigen Test (Lateral Flow Method)	Positive	43	0	43		
	Negative	5	335	340		Negative	2	380	382		
	Total	71	340	411		Total	45	380	425		
Sensitivity: 92.96%		Specificity: 98.53%		Total agreement: 97.57%		Sensitivity: 95.56%		Specificity: 100%		Total agreement: 99.53%	
* NOTE : In field clinical evaluations conducted in Germany, US and UK by lay users.											
CT Value		CT ≤25		CT ≤30		CT ≤33		CT ≤35			
Sensitivity		100%		98.02%		95.88%		92.83%			

W630-Influenza/2019-nCoV Antigen Combo Test	2019-nCoV Antigen	Influenza A	Influenza B
Sensitivity	98.11%	95.81%	94.37%
Specificity	99.72%	91.03%	99.72%
Total Agreement	99.13%	91.86%	99.32%

ORDER INFORMATION

Catalog No.	Product Name	Classification	Packing Size	Sample Type	Storage Condition	Qualification
W196	2019-nCoV Antigen Test (Lateral Flow Method)	Professional Test	20T	Nasopharyngeal swab/oropharyngeal swab	2~30°C	CE
W634	2019-nCoV Antigen Test (Lateral Flow Method)	Self Test	1/5T	Nasal swab	2~30°C	CE 0123
W630	Influenza/2019-nCoV Antigen Combo Test	Professional Test	20T	Nasopharyngeal swab	2~30°C	CE
W771	Flu A&B/2019-nCoV/RSV/ADV/MP Antigen Combo Test	Professional Test	10/20T	Nasopharyngeal swab/oropharyngeal swab	2~30°C	CE
W770	Flu A&B/2019-nCoV/RSV/ADV Antigen Combo Test	Professional Test	10/20T	Nasopharyngeal swab/oropharyngeal swab	2~30°C	CE

Wondfo
Racing For Life

Wondfo

WONDFO
2019-nCoV
ANTIGEN
TEST SOLUTION

RAPID DEFENSE:
CONQUERING COVID-19 AND BEYOND!

- Self test
- Professional test
- Combo test



Guangzhou Wondfo Biotech Co., Ltd.
No.8 Lizhishan Road, Science City,
Huangpu District, 510663, Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com
E-mail: sales@wondfo.com.cn



© 2024 Guangzhou Wondfo Biotech Co., Ltd.
All rights reserved

Guangzhou Wondfo Biotech Co., Ltd.

WONDFO 2019-nCoV ANTIGEN TEST SOLUTION



Direct detection
of the virus



Instant results
within 15mins



Easy to use, no
equipment required



Room temperature
storage(2~30°C)




Sample type: nasopharyngeal,
oropharyngeal or nasal swab



Early detection of COVID-19
(WHO recommends the testing period
is from 3 days before to 5-7 days
after symptoms onset)

CURRENT DIAGNOSTIC METHODS FOR COVID-19



Antigen test

Detect the antigen of the virus, indicating the active viral infection.

RT-PCR

Detect the RNA of virus, indicating the active viral infection.


Antibody test

Detect the antibody generated by immune response after viral infection, indicating the active or past viral infection.

WHEN TO USE ANTIGEN TEST?

Releasing profile

Levels of 2019-nCoV virus and antibodies after infection
*For illustrative purpose only



* Incubation Period: 1~14 days, mostly 3~7 days
* Antibody Window Period: 5-10 days after onset of symptoms

* The minimum detectable limit varies with methodology and sensitivity of test
* IgG antibody test can be referred as one of discharge criteria for recovering COVID-19 patients


ANTIGEN TEST ADVANTAGES

Antigen test OVER RT-PCR

- Short turn-around time (Antigen test: 15 min VS. RT-PCR: 2 hours)
- Inexpensive cost, no equipment required and simple operation make antigen test suitable for point-of-care (POC) setting usage.

Antigen test OVER Antibody test

- Detect the virus directly, allowing the early detection of COVID-19
- Non-invasive sampling (Sampling: swab VS. Blood)




ANTIGEN TEST APPLICATION

Similar to RT-PCR, the detection of antigen indicates the active infection. Under the circumstance that the area(s) is still undergoing widespread community transmission and/or the area with limited RT-PCR resources, antigen can be used for aiding in the diagnosis of COVID-19 suspect patients.


* American CDC also recommends to use rapid antigen tests for screening testing in high-risk congregate settings where the immediate result is required.

Result interpretation



POSITIVE

The patient is undergo active 2019-nCoV infection. Further isolation is required.




NEGATIVE

The patient should be further evaluated by RT-PCR, especially if the result of the antigen test is inconsistent with the clinical context.

ANTIGEN TEST OFFICIAL GUIDELINES



Antigen-detection in the diagnosis of novel coronavirus (2019-nCoV) infection using rapid immunoassays



NPV- negative predictive value PPV- positive predictive value
*The value for NPV and PPV is decided based on products performance and disease prevalence in applied scenarios.

Other antigen test related guidelines

- Interim Guidance for Rapid Antigen Testing for 2019-nCoV, American CDC (8-16-20)
- Considerations for Use of 2019-nCoV Antigen Testing in Nursing Homes, American CDC (8-27-20)
- Antigen-Detection in the Diagnosis of 2019-nCoV Infection Using Rapid Immunoassays, WHO (9-11-20)
- Considerations for Implementation of 2019-nCoV Rapid Antigen Testing, APHL (9-2-20)