

HCG One Step HCG Urine Test

Wondfo One Step HCG Urine Test is a rapid, qualitative detection of human chorionic gonadotropin ("hCG") which appears in urine very early during pregnancy.



RELIABLE

High-quality service is provided.

CONVENIENT


Easy to use.

RAPID

Fast result within 3 minutes.

TEST PROCEDURE & INTERPRETATION


Strip




Cassette



Midstream



ORDER INFORMATION

Catalog No.	Format	Width	Specimen	Sensitivity	Test time	Storage temperature	Shelf life	Qualification
W1-S/ W1-S(20mIU)/ W1-S(10mIU)	Strip	3.5mm	Urine	25/20/10 mIU/mL	3 min	4-30°C	36 months	 0123 FDA (Excluding 20 mIU/mL)
W1-C/ W1-C/ W1-C(20mIU)/ W1-C(10mIU)/ W1-C(10mIU)	Cassette	3.5/4.0mm						
W1-M/ W1-M/ W1-M(20mIU)/ W1-M(10mIU)/ W1-M(10mIU)/ W1-M(10mIU)/ W1-M(10mIU)/ W1-M(10mIU)/ W1-M(10mIU)	Midstream	3.5/5.0/6.0mm						

Guangzhou Wondfo Biotech Co., Ltd.

No.8 Lizhishan Road, Science City, Huangpu District, 510663, Guangzhou, P.R.China

Tel: (+86)400-830-8768

E-mail: sales@wondfo.com.cn

Website: en.wondfo.com

Wondfo reserves all rights. All contents of this brochure are copyrighted by Wondfo.

